

Cy5-PEG-FA (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5-PEG-FA (MW 2000) is a fluorescently labeled biomolecule consisting of a Cy5 dye, polyethylene glycol (PEG), and Folic acid (FA). The Cy5 fluorophore is commonly utilized in applications such as immunolabeling and nucleic acid labeling, with a maximum emission wavelength of approximately 670 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481