


## CRBN ligand-16

## Chemical Properties

CAS No. :	3050687-81-2
Formula:	C <sub>12</sub> H <sub>12</sub> BrNO <sub>2</sub>
Molecular Weight:	282.133
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	CRBN ligand-16 is a CRBN type E3 ubiquitin ligase ligand that can be used in the preparation of PROTAC.
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5445 mL	17.7223 mL	35.4447 mL
5 mM	0.7089 mL	3.5445 mL	7.0889 mL
10 mM	0.3544 mL	1.7722 mL	3.5445 mL
50 mM	0.0709 mL	0.3544 mL	0.7089 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481