

## CRBN ligand-63

## Chemical Properties

CAS No. : 3049446-57-0


Formula: C<sub>12</sub>H<sub>12</sub>BrNO<sub>3</sub>

Molecular Weight: 298.133

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CRBN ligand-63 is a ligand for the cereblon (CRBN) type E3 ubiquitin ligase and can be utilized in the creation of PROTACs.
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3542 mL	16.7712 mL	33.5424 mL
5 mM	0.6708 mL	3.3542 mL	6.7085 mL
10 mM	0.3354 mL	1.6771 mL	3.3542 mL
50 mM	0.0671 mL	0.3354 mL	0.6708 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481