

MRG003

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | MRG003 is an antibody-drug conjugate (ADC) composed of the humanized anti-EGFR IgG1 monoclonal antibody, Becotatug, coupled with MMAE. These components are linked through a valine-citrulline (valine-citrulline) connector, forming the Drug-Linker conjugate VcMMAE within the ADC. |
| Targets(IC50) | EGFR,Antibody-Drug Conjugates (ADCs)   |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481