


Argyrin C

Chemical Properties

CAS No. :	444300-77-0
Formula:	C ₄₁ H ₄₆ N ₁₀ O ₈ S
Molecular Weight:	838.931
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Argyrin C is a natural cyclic nonribosomal peptide that exhibits significant antibacterial activity against <i>Pseudomonas aeruginosa</i> .
Targets(IC50)	Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.192 mL	5.960 mL	11.9199 mL
5 mM	0.2384 mL	1.192 mL	2.384 mL
10 mM	0.1192 mL	0.596 mL	1.192 mL
50 mM	0.0238 mL	0.1192 mL	0.2384 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481