


Tschimganin

Chemical Properties

CAS No. :	38970-49-9
Formula:	C ₁₈ H ₂₄ O ₄
Molecular Weight:	304.381
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tschimganin is a terpene ester with estrogenic activity, found in <i>F. tschimganica</i> . It can stimulate uterine growth in immature rats.
Targets(IC50)	Estrogen Receptor/ERR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2854 mL	16.4268 mL	32.8537 mL
5 mM	0.6571 mL	3.2854 mL	6.5707 mL
10 mM	0.3285 mL	1.6427 mL	3.2854 mL
50 mM	0.0657 mL	0.3285 mL	0.6571 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481