

TRITC-lysine-dextran (MW 70kDa)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TRITC-lysine-dextran, 70kDa, is a dextran conjugated with lysine and contains the fluorescent group TRITC. Lysine-dextran exhibits greater glutamic acid-plasminogen binding capacity compared to dextran alone.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481