

IBI-334

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	IBI-334 is a bispecific antibody targeting both B7-H3 and EGFR. It features an EGFR arm responsible for signal blocking, connected to a modified B7-H3 arm that ensures optimal affinity and binding. The antibody is afucosylated to enhance its antibody-dependent cellular cytotoxicity (ADCC) effects. IBI-334 is widely applicable in EGFR-driven solid tumors.
Targets(IC50)	EGFR,Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481