

Euryachin E

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Euryachin E is an antimicrobial and antiviral agent identified from Eurya chinensis. It shows a MIC of 6.25 µg/mL against B. cereus, a MIC of 0.78 µg/mL against S. aureus, and an IC50 of 6.41 µM against COVID-19. Euryachin E is applicable in research related to infections and inflammation. |
| Targets(IC50) | Antibacterial,Antibiotic,SARS-CoV |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481