


## CRBN ligand-181

## Chemical Properties

CAS No. :	3036607-12-9
Formula:	C <sub>12</sub> H <sub>11</sub> BrFNO <sub>2</sub>
Molecular Weight:	300.124
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	CRBN ligand-181 is a CRBN-type E3 ubiquitin ligase ligand, suitable for use in PROTAC preparation.
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.332 mL	16.660 mL	33.320 mL
5 mM	0.6664 mL	3.332 mL	6.664 mL
10 mM	0.3332 mL	1.666 mL	3.332 mL
50 mM	0.0666 mL	0.3332 mL	0.6664 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481