

## Red Fluorescent PLGA nanoparticles, 500nm

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Red Fluorescent PLGA nanoparticles, 500nm, are nanoparticles composed of the polymer PLGA and tagged with a red fluorescent dye. These nanoparticles are suitable for applications in drug delivery and tracing.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481