


## Lenalidomide-hex-5-ynoic acid

## Chemical Properties

CAS No. :	2379870-47-8
Formula:	C <sub>19</sub> H <sub>18</sub> N <sub>2</sub> O <sub>5</sub>
Molecular Weight:	354.357
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Lenalidomide-hex-5-ynoic acid is a conjugate of an E3 ligase ligand and linker, utilized in the synthesis of PROTACs, such as PROTACcGAS degrader-1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.822 mL	14.1099 mL	28.2199 mL
5 mM	0.5644 mL	2.822 mL	5.644 mL
10 mM	0.2822 mL	1.411 mL	2.822 mL
50 mM	0.0564 mL	0.2822 mL	0.5644 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481