


## LC3B ligand 1

## Chemical Properties

CAS No. :	902959-52-8
Formula:	C <sub>14</sub> H <sub>11</sub> N <sub>3</sub> O <sub>2</sub>
Molecular Weight:	253.256
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	LC3B ligand 1 is the target protein ligand for PROTAC, and it can be utilized in the synthesis of PROTACCDK9 autophagic degrader 1.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9485 mL	19.7426 mL	39.4851 mL
5 mM	0.7897 mL	3.9485 mL	7.897 mL
10 mM	0.3949 mL	1.9743 mL	3.9485 mL
50 mM	0.079 mL	0.3949 mL	0.7897 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481