

MAS825

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	MAS825 is a bispecific monoclonal antibody targeting IL-1 β and IL-18, designed to reduce inflammasome-mediated inflammation. This compound combines the inhibitory domains of anti-IL-1 β and anti-IL-18 monoclonal antibodies (mAb) into one high-affinity cytokine-capturing bispecific molecule.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481