


## Aberrant tau degrader 2

## Chemical Properties

CAS No. :	2417159-58-9
Formula:	C19H27N7O3S
Molecular Weight:	433.528
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Aberrant tau degrader 2 (Compound 4-12) acts as a degrader of tau protein and serves as a target protein ligand for the synthesis of PROTAC degrader C004019.
Targets(IC50)	Microtubule Associated,Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3066 mL	11.5332 mL	23.0665 mL
5 mM	0.4613 mL	2.3066 mL	4.6133 mL
10 mM	0.2307 mL	1.1533 mL	2.3066 mL
50 mM	0.0461 mL	0.2307 mL	0.4613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481