

## Anti-TrkB/NTRK2 Antibody

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-TrkB/NTRK2 Antibody is a highly specialized monoclonal antibody target the extracellular domain of the tropomyosin receptor kinase B (TrkB/NTRK2), modulation of neurotrophic signaling pathways to promote neuronal survival and synaptic plasticity across various preclinical neurological experimental models.
Targets(IC50)	Trk receptor
In vitro	In reporter cell line assays, Anti-TrkB/NTRK2 Antibody (1C2, 38B8) induced receptor phosphorylation and activated intracellular signal transduction pathways [1].

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481