


## WDR5 ligand 2

## Chemical Properties

CAS No. :	2737222-78-3
Formula:	C <sub>29</sub> H <sub>31</sub> F <sub>3</sub> N <sub>4</sub> O <sub>4</sub>
Molecular Weight:	556.576
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	WDR5ligand 2 is a ligand for WDR5 and can be utilized in the synthesis of PROTAC WDR5degrader 1.
Targets(IC50)	JAK,Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7967 mL	8.9834 mL	17.9669 mL
5 mM	0.3593 mL	1.7967 mL	3.5934 mL
10 mM	0.1797 mL	0.8983 mL	1.7967 mL
50 mM	0.0359 mL	0.1797 mL	0.3593 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481