


Atrovenetin

Chemical Properties

CAS No. :	2582-86-7
Formula:	C ₁₉ H ₁₈ O ₆
Molecular Weight:	342.343
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Atrovenetin is an antibiotic found in <i>Penicillium</i> species (antibiotic). It exhibits inhibitory activity against <i>Bacillus subtilis</i> and <i>Staphylococcus aureus</i> . Atrovenetin is also a potent antioxidant.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9211 mL	14.6054 mL	29.2107 mL
5 mM	0.5842 mL	2.9211 mL	5.8421 mL
10 mM	0.2921 mL	1.4605 mL	2.9211 mL
50 mM	0.0584 mL	0.2921 mL	0.5842 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481