

CypA ligand-2

Chemical Properties

CAS No. :

Formula: C₃₄H₃₄N₆O₇

Molecular Weight: 638.67

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CypA ligand-2 serves as the target protein ligand for PROTAC RJS308. RJS308 is a PROTAC specifically designed to target CypA, demonstrating antiviral properties.
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5658 mL	7.8288 mL	15.6575 mL
5 mM	0.3132 mL	1.5658 mL	3.1315 mL
10 mM	0.1566 mL	0.7829 mL	1.5658 mL
50 mM	0.0313 mL	0.1566 mL	0.3132 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481