

Cy5.5-PEG-SH (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5.5-PEG-SH (MW 3400) is a fluorescent labeling reagent consisting of the Cy5.5 fluorescent dye, polyethylene glycol (PEG), and a thiol (SH) group, making it suitable for molecular conjugation. The SH group is highly reactive, allowing it to form stable thioether bonds with various functional groups, such as maleimide.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481