

Cy5.5-PEG-SH (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5.5-PEG-SH (MW 10000) is a fluorescent labeling reagent that includes Cy5.5 fluorescent dye, polyethylene glycol (PEG), and a thiol group (SH), making it suitable for molecular conjugation. The highly reactive SH group can form stable thioether bonds by reacting with various functional groups such as maleimide.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481