

IMC-D11

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	IMC-D11 is a fully human monoclonal antibody that specifically targets the extracellular domain of fibroblast growth factor receptor 3 (FGFR3) and blocks ligand-induced activation. By inhibiting FGFR3 phosphorylation and its downstream signaling cascades, this antibody effectively suppresses the growth and survival of various FGFR3-dependent tumor cells.
Targets(IC50)	FGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481