


## E3 ligase Ligand 56

## Chemical Properties

CAS No. :	1472731-41-1
Formula:	C11H13N3O2
Molecular Weight:	219.24
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	E3 Ligase Ligand 56 is a ligand for the E3 ligase Cereblon, and it can be utilized in the synthesis of [PROTACCDK9degrader-11].
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5612 mL	22.8061 mL	45.6121 mL
5 mM	0.9122 mL	4.5612 mL	9.1224 mL
10 mM	0.4561 mL	2.2806 mL	4.5612 mL
50 mM	0.0912 mL	0.4561 mL	0.9122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481