

Cy5.5-PEG-Biotin (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5.5-PEG-Biotin (MW 5000) is a polyethylene glycol derivative that incorporates Cy5.5 fluorescent dye, polyethylene glycol (PEG), and biotin. The Cy5.5 fluorophore is commonly utilized in applications such as immunolabeling and nucleic acid labeling. Biotin exhibits a high affinity for binding with streptavidin or avidin. In the specific capture and detection of biomolecules, Cy5.5-PEG-Biotin (MW 5000) plays a crucial role.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481