


Lomonitinib

Chemical Properties

CAS No. :	2923221-56-9
Formula:	C ₂₇ H ₂₄ N ₄ O ₂
Molecular Weight:	436.505
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lomonitinib is a potent and selective pan-FLT3/IRAK4 inhibitor with antitumor properties. It shows promise for research in myeloid leukemia.
Targets(IC50)	FLT,IRAK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2909 mL	11.4545 mL	22.909 mL
5 mM	0.4582 mL	2.2909 mL	4.5818 mL
10 mM	0.2291 mL	1.1454 mL	2.2909 mL
50 mM	0.0458 mL	0.2291 mL	0.4582 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481