


CRBN ligand-83

Chemical Properties

CAS No. :	3036607-28-7
Formula:	C ₁₂ H ₁₁ BrFNO ₃
Molecular Weight:	316.123
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	CRBN ligand-83 is a ligand of the CRBN-type E3 ubiquitin ligase, suitable for the preparation of PROTACs.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1634 mL	15.8168 mL	31.6336 mL
5 mM	0.6327 mL	3.1634 mL	6.3267 mL
10 mM	0.3163 mL	1.5817 mL	3.1634 mL
50 mM	0.0633 mL	0.3163 mL	0.6327 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481