

PM-8002

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PM-8002 is a bispecific antibody that targets PD-L1 and VEGF-A. It is applicable for research on solid tumors.
Targets(IC50)	PD-1/PD-L1,VEGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481