


p38 Kinase inhibitor 8

Chemical Properties

CAS No. :	815595-25-6
Formula:	C ₂₂ H ₂₁ FN ₆ O ₂
Molecular Weight:	420.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	p38 Kinase inhibitor 8 (Compound CCLXXVIII) is an orally active inhibitor targeting p38 β and JNK2 α 2, with IC ₅₀ values of 6.3 nM and 53.6 nM, respectively. It has demonstrated anti-inflammatory effects in a rat model of collagen-induced arthritis.
Targets(IC ₅₀)	JNK,p38 MAPK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3785 mL	11.8923 mL	23.7846 mL
5 mM	0.4757 mL	2.3785 mL	4.7569 mL
10 mM	0.2378 mL	1.1892 mL	2.3785 mL
50 mM	0.0476 mL	0.2378 mL	0.4757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481