

MGD010

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	MGD010 is a bispecific antibody targeting CD32B (FcγRIIb) and CD79B. By simultaneously co-engaging CD32B (FcγRIIb) and the BCR component CD79B, MGD010 delivers an inhibitory signal, thereby suppressing activated B cells. MGD010 shows potential for research in autoimmune diseases.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481