


Argyrin G

Chemical Properties

CAS No. :	444300-75-8
Formula:	C41H46N10O9S
Molecular Weight:	854.931
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Argyrin G is a naturally occurring cyclic nonribosomal peptide.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1697 mL	5.8484 mL	11.6969 mL
5 mM	0.2339 mL	1.1697 mL	2.3394 mL
10 mM	0.117 mL	0.5848 mL	1.1697 mL
50 mM	0.0234 mL	0.117 mL	0.2339 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481