

DT-3 acetate

Chemical Properties

CAS No. :

Formula: C152H258N52O28S.xC2H4O2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DT-3 acetate is a cell-permeable peptide that acts as an inhibitor of protein kinase G (PKG α), effectively blocking the cGMP-PKG signaling pathway.
Targets(IC50)	PKA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481