

TRITC-lysine-dextran (MW 10kDa)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	TRITC-lysine-dextran, 10kDa, is a dextran compound conjugated with lysine and contains the fluorescent group TRITC. Lysine-dextran exhibits a higher binding affinity to glutamic acid-plasminogen compared to dextran alone.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481