

ZSNI-21

Chemical Properties

CAS No. :

Formula: C₂₆H₂₅N₃O₅

Molecular Weight: 459.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ZSNI-21 is a dual inhibitor targeting ADAM17 and HDAC2. It effectively suppresses the proliferation of Bel-7402 cells and exhibits significant anti-metastatic properties against HCC-LM3 cells, making it a promising candidate for hepatocellular carcinoma (HCC) research.
Targets(IC50)	HDAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1763 mL	10.8816 mL	21.7633 mL
5 mM	0.4353 mL	2.1763 mL	4.3527 mL
10 mM	0.2176 mL	1.0882 mL	2.1763 mL
50 mM	0.0435 mL	0.2176 mL	0.4353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481