

## Haptoglobin, Phenotype 2-2, Human Plasma

### Chemical Properties

CAS No. : 9087-69-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Haptoglobin, Phenotype 2-2 is a plasma protein that binds hemoglobin (Hb) to remove free hemoglobin and possesses antioxidant activity. It can be utilized in research related to kidney failure and cardiovascular diseases associated with diabetes.
Targets(IC50)	Endogenous Metabolite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481