

1,2-Isopropylidene-3-oleoyl-sn-glycerol

Chemical Properties


CAS No. : 33001-45-5

Formula: C₂₄H₄₄O₄

Molecular Weight: 396.604

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,2-Isopropylidene-3-oleoyl-sn-glycerol is a fatty acid that serves as a biomaterial or organic compound utilized in life science research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5214 mL	12.6072 mL	25.2143 mL
5 mM	0.5043 mL	2.5214 mL	5.0429 mL
10 mM	0.2521 mL	1.2607 mL	2.5214 mL
50 mM	0.0504 mL	0.2521 mL	0.5043 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481