

Cy5-PEG-biotin (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5-PEG-biotin (MW 5000) is a polyethylene glycol derivative that incorporates Cy5 fluorescent dye, polyethylene glycol (PEG), and biotin. The Cy5 fluorophore is commonly utilized for applications such as immunolabeling and nucleic acid labeling. Biotin offers high-affinity binding with streptavidin or avidin. Cy5-PEG-biotin (MW 5000) plays a crucial role in the specific capture and detection of biomolecules.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481