


## 3-(4-Bromo-2-chlorophenyl)piperidine-2,6-dione

## Chemical Properties

CAS No. :	2923549-60-2
Formula:	C <sub>11</sub> H <sub>9</sub> BrClNO <sub>2</sub>
Molecular Weight:	302.552
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	3-(4-Bromo-2-chlorophenyl)piperidine-2,6-dione acts as a ligand for CRBN type E3 ubiquitin ligase and is useful in the synthesis of PROTAC.
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3052 mL	16.5262 mL	33.0524 mL
5 mM	0.661 mL	3.3052 mL	6.6105 mL
10 mM	0.3305 mL	1.6526 mL	3.3052 mL
50 mM	0.0661 mL	0.3305 mL	0.661 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481