


Fenfuram

Chemical Properties

CAS No. :	24691-80-3
Formula:	C ₁₂ H ₁₁ NO ₂
Molecular Weight:	201.221
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fenfuram is an antifungal agent with an EC ₅₀ of 6.18 mg/mL against R. solani. It is utilized in the research of anti-infective treatments.
Targets(IC ₅₀)	Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9697 mL	24.8484 mL	49.6968 mL
5 mM	0.9939 mL	4.9697 mL	9.9394 mL
10 mM	0.497 mL	2.4848 mL	4.9697 mL
50 mM	0.0994 mL	0.497 mL	0.9939 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481