


Asperaldin

Chemical Properties

CAS No. : 561297-46-9
 Formula: C₁₆H₁₈O₅
 Molecular Weight: 290.311
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Asperaldin is a potent inhibitor of aldose reductase (RLAR) and can be found in <i>Aspergillus niger</i> CFR-1046.
Targets(IC50)	Reductase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4446 mL	17.223 mL	34.4459 mL
5 mM	0.6889 mL	3.4446 mL	6.8892 mL
10 mM	0.3445 mL	1.7223 mL	3.4446 mL
50 mM	0.0689 mL	0.3445 mL	0.6889 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481