

Anti-human IL-25 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-human IL-25 Antibody is a specific monoclonal antibody directed against the pro-inflammatory cytokine interleukin-25 (IL-25, also known as IL-17E). By neutralizing circulating IL-25, this antibody blocks the overactivation of Th2-type immune responses, playing a critical role in the study of allergic asthma, dermatitis, and related inflammatory conditions.
Targets(IC50)	IL Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481