


E3 Ligase Ligand-linker Conjugate 165

Chemical Properties

CAS No. :	2789679-85-0
Formula:	C ₂₁ H ₃₀ N ₄ O ₂
Molecular Weight:	370.488
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	E3 Ligase Ligand-linker Conjugate 165 is a conjugate of an E3 ligase ligand and a linker, suitable for synthesizing PROTACS J44236.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6991 mL	13.4956 mL	26.9913 mL
5 mM	0.5398 mL	2.6991 mL	5.3983 mL
10 mM	0.2699 mL	1.3496 mL	2.6991 mL
50 mM	0.054 mL	0.2699 mL	0.5398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481