


## CRBN ligand-74

## Chemical Properties

CAS No. :	3065251-87-5
Formula:	C <sub>11</sub> H <sub>8</sub> BrF <sub>2</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	304.088
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	CRBN ligand-74 is a CRBN-type E3 ubiquitin ligase ligand useful for the preparation of PROTACs (Proteolysis Targeting Chimeras).
Targets(IC50)	Ligands for E3 Ligase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2885 mL	16.4425 mL	32.885 mL
5 mM	0.6577 mL	3.2885 mL	6.577 mL
10 mM	0.3289 mL	1.6443 mL	3.2885 mL
50 mM	0.0658 mL	0.3289 mL	0.6577 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481