

Lapemelanotide zapixetan

Chemical Properties

CAS No. :

Formula: C79H119N25O18

Molecular Weight: 1706.95

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Lapemelanotide zapixetan exhibits antitumor activity.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5858 mL	2.9292 mL	5.8584 mL
5 mM	0.1172 mL	0.5858 mL	1.1717 mL
10 mM	0.0586 mL	0.2929 mL	0.5858 mL
50 mM	0.0117 mL	0.0586 mL	0.1172 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481