


FGFR1 inhibitor-15

Chemical Properties

CAS No. :	459137-97-4
Formula:	C17H13FN4O
Molecular Weight:	308.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	FGFR1 inhibitor-15 (Compound 23) is an FGFR1 inhibitor with an IC50 value of 27 μ M, useful for tumor research.
Targets(IC50)	FGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2435 mL	16.2174 mL	32.4349 mL
5 mM	0.6487 mL	3.2435 mL	6.487 mL
10 mM	0.3243 mL	1.6217 mL	3.2435 mL
50 mM	0.0649 mL	0.3243 mL	0.6487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481