

Bovine Protein C

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bovine Protein C is activated into the serine protease activated protein C (APC) by α -thrombin or the α -thrombin/thrombomodulin complex. It selectively inactivates factors Va and VIIIa, functioning as a potent anticoagulant.
Targets(IC50)	APC/C

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481