


GW842166X

Chemical Properties

CAS No. : 666260-75-9
 Formula: C₁₈H₁₇Cl₂F₃N₄O₂
 Molecular Weight: 449.25
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	GW842166X is a potent and highly selective agonist of cannabinoid receptor CB2 receptor with EC ₅₀ of 63 nM, shows no significant activity at CB1 receptor. Phase 2.
Targets(IC ₅₀)	Cannabinoid Receptor
In vitro	GW842166X demonstrates comparable potency and efficacy for rat and human recombinant CB2 receptors, with EC ₅₀ values of 91 nM and 63 nM, respectively. [1] It exhibits full agonist activity in cyclase assays with an EC ₅₀ of 133 nM and an E _{max} of 101%, and weak agonist activity in FLIPR assays with an EC ₅₀ of 7.780 μM and an E _{max} of 84%. [2]
In vivo	GW842166X has an oral bioavailability of 58% and a half-life of 3 h when dosed orally in the rat. GW842166X has extremely high potency with an oral ED ₅₀ of 0.1 mg/kg and shows full reversal of hyperalgesia at 0.3 mg/kg in the FCAa model of inflammatory pain. [1] GW842166X orally administrated at a dose of 15 mg/kg for 8 days produced a significant reversal of the CCI induced decrease in paw withdrawal threshold in a rat model of neuropathic pain. [3]

Solubility Information

Solubility	H ₂ O: < 1 mg/mL (insoluble or slightly soluble), Ethanol: < 1 mg/mL (insoluble or slightly soluble), DMSO: 16 mg/mL (35.61 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2259 mL	11.1297 mL	22.2593 mL
5 mM	0.4452 mL	2.2259 mL	4.4519 mL
10 mM	0.2226 mL	1.113 mL	2.2259 mL
50 mM	0.0445 mL	0.2226 mL	0.4452 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Giblin GM, et al. J Med Chem, 2007, 50(11), 2597-2600.

He Y, Jia H, Yang Q, et al. Specific Activation of CB2R Ameliorates Psoriasis-Like Skin Lesions by Inhibiting Inflammation and Oxidative Stress. Inflammation. 2023: 1-17.

Yao BB, et al. Br J Pharmacol, 2008, 153(2), 390-401.

Clayton NM, et al. Proceedings of the British Pharmacological Society at <https://www.pA2online.org/Vol2Issue4abst050P.html>

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481