

SGC 6870N

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Negative control for SGC 6870. SGC 6870N is the (S)-enantiomer, which is inactive on PRMT6.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 46.94 mg/mL, Sonication is recommended. Ethanol: 46.94 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481