

## Magainin Spacer Peptide

## Chemical Properties

CAS No. :


Formula: C97H150N28O39

Molecular Weight: 2332.39

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Magainin Spacer Peptide is a polypeptide molecule with the sequence DAEAVGPEAFADQDLDEREVR.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4287 mL	2.1437 mL	4.2874 mL
5 mM	0.0857 mL	0.4287 mL	0.8575 mL
10 mM	0.0429 mL	0.2144 mL	0.4287 mL
50 mM	0.0086 mL	0.0429 mL	0.0857 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481