


## [Tyr0] Bradykinin

## Chemical Properties

CAS No. :	33289-76-8
Formula:	C59H82N16O13
Molecular Weight:	1223.38
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	[Tyr0] Bradykinin is a polypeptide molecule with the sequence YRPPGFSPFR.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8174 mL	4.087 mL	8.1741 mL
5 mM	0.1635 mL	0.8174 mL	1.6348 mL
10 mM	0.0817 mL	0.4087 mL	0.8174 mL
50 mM	0.0163 mL	0.0817 mL	0.1635 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481