


## Fam-Bradykinin

## Chemical Properties

CAS No. :	5979-11-3
Formula:	C <sub>56</sub> H <sub>85</sub> N <sub>15</sub> O <sub>17</sub>
Molecular Weight:	1240.36
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Fam-Bradykinin is a polypeptide molecule with the sequence 5FAM-RPPGFSPFR.
Targets(IC50)	Endogenous Metabolite, Bradykinin Receptor, Serine Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8062 mL	4.0311 mL	8.0622 mL
5 mM	0.1612 mL	0.8062 mL	1.6124 mL
10 mM	0.0806 mL	0.4031 mL	0.8062 mL
50 mM	0.0161 mL	0.0806 mL	0.1612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481